



MAGic screen magnification software is an innovative low vision solution for computer users. By combining state-of-the-art magnification features with a screen reading speech option, MAGic helps people with low vision view information on their computer screen while hearing it through a built-in speech synthesizer. Use magnification and speech together or independently.

## ***MAGic at a Glance***

- Magnified on-screen images combined with a speech option provide easier access to computer applications
- 45 magnification levels deliver magnification of from 1x to 36x original size
- Hear on-screen text as it is read to you with powerful, responsive, and customizable low vision speech
- Use screen magnification and the speech option together or independently
- Tailor the on-screen appearance in ways that best address your needs and preferences

## **Customizable Views**

With 45 magnification levels, MAGic magnifies the contents of your computer screen up to 36 times its actual size. Powerful color enhancement features give you a choice of many low-glare, high-contrast combinations. Apply tinting, view the screen as two colors (monochrome), invert brightness and colors, and switch colors. These enhancements can be applied individually or in combination to precisely tailor the screen's appearance.

Use full-screen mode or see both magnified and unmagnified areas of your screen simultaneously with partial-screen magnification modes: Lens View, Overlay View, Split View, and Dynamic Lens View.

## Tracking Enhancements

With the unique tracking enhancement features of MAGic, follow the movement of your mouse, your cursor, and the application focus on the screen. Change the color of these enhancements, apply transparency effects, adjust their size, and more to make navigating screens quick and easy.

## Speech Features

MAGic has the flexibility of two included text-to-speech software synthesizers: RealSpeak™ Solo and Eloquence. Choose the one that is right for your use. MAGic also provides a variety of speech settings that allow you to control the amount of speech you hear and the properties of the voices. Whether navigating by keyboard or mouse, you'll get quick access to information as MAGic announces application events and follows your navigation of the screen. Easily clarify words or browse a document by moving the mouse or cursor to the text you want spoken. Navigate and read by word, line, sentence, or paragraph using the built-in functions.

## Enhanced Productivity

The smooth, seamless performance of MAGic and its array of flexible options help you get the most from your computer. You maximize your efficiency and productivity because MAGic automatically tracks the activities of your applications.

More than five hours of built-in training tools make it easy to learn MAGic,

including audio training materials, talking installation, easy-to-use online Help, and context-sensitive Help for all MAGic controls. The MAGic Context Menu allows you to change preferences instantly without having to remember keystrokes or return to the main window.

## Additional Key Features Include:

- Remote access via Citrix and Terminal Services (extra authorization is required)
- Available multi-user license agreements
- Specific settings created by users are automatically applied when you switch between applications
- Full compatibility and handshaking with JAWS® for Windows®, the world's most popular screen reader, enables you to combine the powerful features of MAGic and JAWS to create a dynamic vision solution
- MAGic EyeMerge™, an image management feature that provides seamless integration with our ONYX™ multiple-view cameras, is included

## Contact your local dealer or contact us:

Freedom Scientific  
Blind/Low Vision Group  
11800 31<sup>st</sup> Court North  
St. Petersburg, Florida 33716-1805  
Phone/Sales: 800-444-4443  
727-803-8000

[www.FreedomScientific.com](http://www.FreedomScientific.com)

