

Eye-Pal[®] Ace Plus

Portable Scanning and Reading Solution



Accurately and Instantly Read and Magnify Text!

Easy-to-use, portable, and lightweight - for the home, office, and on the go!



Features

- ▶ Read content aloud
- ▶ Built-in WiFi and email
- ▶ Set reminders with AudioMinder
- ▶ Download books
- ▶ Compatible with Focus Braille Displays

Learn more about getting back to the things that matter most.

www.FreedomScientific.com



WiFi enabled



Exclusive AudioMinder technology



1-866-569-8474 • 727-803-8000 (outside US)

Eye-Pal[®] Ace Plus

Portable Scanning and Reading Solution



Remain Active and Stay Connected!

Browse the web, check email, and clearly view menus, prescriptions, and more!

Built-in thumbwheels, large buttons, and anti-glare lights make operation easy, even in low-light situations. With less than 3 seconds from scanning to reading, the Eye-Pal Ace Plus will keep up with your lifestyle.

Specifications

Display	10-inch (diagonal) color display
Weight	3.5 pounds
Power	Rechargeable battery and/or AC
Controls	Built-in thumbwheels and large buttons
Dimensions	12.4" W x 10.6" H x 3.4" D
Camera	Two 5-MP auto-focus HD cameras
Light source	Dual anti-glare built-in ultra bright LEDs
Energy-saving features	Auto-standby, auto-sleep modes
Language options	20 languages/voices



www.FreedomScientific.com

For product details, including specifications and warranty, please visit our website or scan this code.

750557-001 Rev. A



1-866-569-8474 • 727-803-8000 (outside US)