

ONYX[™]

Freedom Scientific

Deskset XL



ONYX Deskset XL at a Glance:

- Self-contained and easy to carry
- Magnification up to 64x (17-inch monitor) and 71x (19-inch monitor)
- Distance, document, and self views
- 350 degree camera rotation
- Adjustable monitor height and tilt



Power and Productivity

The ONYX Deskset XL is an ideal "on-the-go" low vision solution for clear, powerful magnification that is portable, adaptable, and easy to use. It fits in any classroom, business, or home environment and moves from one viewing situation to the next as quickly as you do.

The telescoping camera arm of the ONYX Deskset XL provides working space under the camera and a vantage point for distance viewing. It collapses for compact transport. The ONYX Deskset XL features the versatile ONYX camera mounted on a rotating arm to a lightweight flat panel monitor. It provides you with the ultimate in distance and up-close magnification in a self-contained unit that plugs into any standard wall socket.

Clear Images – Near and Far

The ONYX Deskset XL produces precise magnification in three distinct viewing modes. In the distance mode, even objects from across large rooms or auditoriums can be seen in sharp focus. The document mode magnifies items such as reading materials, diagrams, product labels, and craft projects. With the self view, you get a magnified true mirror image close up. The camera remembers its settings in each position, so you can switch back and forth between viewing modes without refocusing or readjusting the camera.

Versatility for any Situation

The camera's mounting arm has three axes, and the camera itself rotates





up to 350 degrees, enabling you to place it in an almost endless number of positions. The monitor can be raised and lowered and tilted forward and backward for maximum comfort.

The many key features of the ONYX Deskset XL include:

- A choice of full color, enhanced black on white, enhanced white on black, yellow on blue, and yellow on black. These modes can be used to enhance contrast for reading text
- A full range of brightness adjustments to eliminate glare and enable viewing in low light
- Masking and Reading Lines
- The Find Feature easily establishes viewing perspective on an image prior to zooming in
- Focus Lock keeps an object or document crisp and clear when working or writing under the camera
- Freeze Frame provides ample time to carefully review images
- A unique, pocket-sized RF remote control



Grasp the handle, and slip your ONYX Deskset XL into its rolling case.

Specifications

Weight:

- ONYX Deskset XL 17: 12 lbs / 5.44 kg
- ONYX Deskset XL 19: 15.4 lbs / 7 kg
- Camera-only weight: 14.8 oz / 420 g

Working space

under camera:
17 3/4 in / 450.85 mm

Magnification at

- 17 3/4 in / 450.85 mm:
- 17" monitor: 64x
 - 19" monitor: 71x

Contact your local dealer or contact us:

Freedom Scientific
11800 31st Court North
St. Petersburg, Florida 33716-1805
Phone/Sales: 800-444-4443
www.FreedomScientific.com

About Freedom Scientific

Freedom Scientific is the world's leading manufacturer of assistive technology products for those who are blind or vision impaired. For more information, e-mail info@freedomscientific.com