

RUBY®

Handheld Video Magnifier

RUBY at a Glance

- 4.3-inch high resolution LCD display
- Lightweight design — 7.7 oz (218g)
- 5 viewing modes
- Rapid reading without blur
- Freeze Frame and change magnification
- Save up to 15 screen shots
- Folding handle
- 2 hours continuous use
- Use rechargeable or disposable batteries
- Onscreen battery status indicator on start up



For the many people who need help reading text and seeing details because of low vision, traditional magnifying glasses work only up to a point. The RUBY takes handheld magnification to the next level and beyond. The RUBY's full color, high brightness video screen makes it outstanding for magnifying bills, mail, receipts, and checks. It easily slips into a pocket or purse as the perfect traveling companion for visiting the grocery store, pharmacy, bank, library, restaurant, or anywhere else. And it's quick and easy to operate, with bright, colorful push buttons.

Lay the RUBY flat on a surface, or unfold the handle and use it like a standard magnifying glass — but see everything with enhanced contrast and NO DISTORTION.



 **Freedom
Scientific®**
everything in sight™

Change the text and background color combinations to suit your most comfortable reading. The RUBY features full color mode and four enhanced contrast modes for reading.

Read as fast as you want. The screen image stays crisp and clear, even when rapidly moving the RUBY across lines of text — no motion blur.

Click the Freeze Frame button to capture an image and bring it close for comfortable reading — even change the magnification of the still image. Store up to 15 images, making the RUBY an ideal low vision aid for referring to prescription instructions, phone numbers, and other frequently needed information.





The RUBY comes with a 30-day money-back guarantee.

Magnification Range: 2x to 14x with handle extended; 5x, 7.5x, and 10x flat

Display: 4.3" diagonal (10.9 cm); 16:9 aspect ratio

Contrast Modes

Full color mode and four enhanced contrast modes for reading:

-  Black text on a white background
-  White text on a black background
-  Yellow text on a blue background
-  Yellow text on a black background

Dimensions: 4.9" x 2.9" x 1.35" (.65" thick with handle extended)
12.5 cm x 7.5 cm x 3.4 cm (1.7 cm thick with handle extended)

Weight: 7.7 ounces (218g)

Battery: 4 AAA rechargeable batteries with 2 hours continuous use; 2 hour charge time. Onscreen battery status indicator on startup and flashing LED low battery indicator.

Also works with disposable batteries

Carrying case and charger included

Specifications subject to change

About Freedom Scientific

Freedom Scientific is the world's leading manufacturer of assistive technology products for those who are blind or vision impaired. For more information, e-mail info@FreedomScientific.com



www.FreedomScientific.com