

RUBY™



RUBY at a Glance

- 4.3-inch high resolution LCD display
- Lightweight design - 7.7 oz (218g)
- 5 viewing modes
- Freeze frame and change magnification
- Folding handle
- 2 hours continuous use
- Use rechargeable or disposable batteries

RUBY™ Handheld Video Magnifier

For the many people who need help reading text and seeing details because of low vision, traditional magnifying glasses work only up to a point. The RUBY™ video magnifier takes handheld magnification to the next level and beyond. The RUBY's 4.3-inch full color, high brightness video screen makes it outstanding for reading bills, mail, receipts, and checks. It is so small and unobtrusive that it easily slips into a pocket or purse as the perfect traveling companion for visiting the grocery store, the pharmacy, the bank, the library, bookstore, restaurant, or anywhere else.

The RUBY is easy to use for those who are unfamiliar with technology products. Just turn it on with one button, place it over an object, and adjust enlargement to your preference with the easy zoom button. The unit stays flat on the table, and the magnified image appears in the display window. With only one button, you can adjust the screen from full color for viewing pictures to any of four high contrast text modes.

 **Freedom
Scientific**
everything in **sight**™

www.LowVisionSolutions.com

RUBY™



Handheld Video Magnifier

Flip out the convenient handle and use the RUBY like a conventional magnifying glass – but see everything with enhanced contrast and NO DISTORTION.

Press the Freeze Frame button to capture the image and bring it close for comfortable reading - even change the magnification of the still image.



SPECIFICATIONS

Magnification Range:

2x to 14x with handle extended
5x, 7.5x, and 10x flat on the paper

Camera:

Full color mode and four enhanced contrast modes for reading: black text on a white background, white text on a black background, yellow text on a blue background, or yellow text on a black background

Battery:

4 AAA rechargeable batteries
2 hours continuous use
2 hour charge time
Also works with disposable batteries

Display:

4.3" diagonal (102mm); 16:9 aspect ratio

Dimensions:

4.94" x 2.94" x 1.35" (.65" thick with handle extended)
125 mm x 75 mm x 34 mm (16.5 mm thick with handle extended)

Weight:

7.7 ounces (218g)

Carrying case and charger with universal adapters included

Specifications subject to change

About Freedom Scientific

Freedom Scientific is the world's leading manufacturer of assistive technology products for those who are blind or vision impaired. For more information, e-mail Info@FreedomScientific.com



www.LowVisionSolutions.com