



# NEW PRODUCT!

## RUBY<sup>®</sup> 7 HD

Handheld Video Magnifier



...the newest addition to the RUBY<sup>®</sup> family.

# RUBY<sup>®</sup> 7 HD

Handheld Video Magnifier



**NEW  
PRODUCT!**



## Features

- 2x - 24x magnification
- Large, 7-inch widescreen LCD
- Built-in reading stand
- Lightweight - 18 oz/510 g
- Full size HDMI port - connect to TV
- 3 hours continuous use
- 20 high-contrast colors
- Freeze Frame
- Reading Line and Masks
- Save up to 100 images
- PivotCam™ Rotating Camera
- Panning button for still images

## What's Included

- Carrying case
- Charger
- USB cable
- HDMI cable
- User guide

## Dimensions

8.3 in x 5.8 in. x 1.0 in  
21.1 cm x 14.7 cm x 2.5 cm



750579-001 Rev. A