

TOPAZ[®] XL

Desktop Video Magnifier with PC Connectivity

Features at a Glance

- Clear, vivid images
- Magnifies up to 77x
- USB and VGA connectivity to your PC
- Freeze Frame, Focus Lock, Find feature, and Position Locator
- Directed LED lighting for a bright work surface without glare
- 8.25 inches of workspace under the camera
- 30 high-contrast text color modes plus true color and enhanced full color
- Monitor raises, lowers, tilts, and swivels
- Unique groove on XY table holds bottles and cans still for viewing



Space-Saving, High-Powered Video Magnifier with PC Connectivity

The TOPAZ XL has all the great ease-of-use features of our other TOPAZ video magnifiers plus two ways to connect to your computer.

Connect the TOPAZ XL to your computer with USB, and use Freedom Scientific's optional GEM™ software to work with your PC and document simultaneously. Split your computer screen vertically to enter a column of numbers in a spreadsheet or horizontally to type notes from your document in a word processor. Even save images from your TOPAZ XL.

Connect your computer's monitor cable to TOPAZ XL's VGA port and use it to view documents as well as your PC. Toggle back and forth with a button push. Save space on your desk, and avoid moving your chair back and forth.

Add Reading Lines and Masks – and More

With built-in masks and reading lines, you can block all but the line you are reading to reduce eye fatigue and set reference lines as you move a document under the camera. You can choose horizontal or vertical markers, move them anywhere on the screen, and adjust the space between.



Learn more online at www.FreedomScientific.com

TOPAZ[®] XL

Desktop Video Magnifier with PC Connectivity

Magnification:

- 17-inch monitor - 1.9x - 49x
- 20-inch monitor - 2.4x - 64x
- 22-inch monitor - 2.6x - 68x
- 24-inch monitor - 2.9x - 77x

Color Modes:

Full Color, Black on White, White on Black, Yellow on Blue, Yellow on Black, plus Enhanced Full Color and 27 additional, configurable color modes

Working Space Under Camera:

8.25 in. (21 cm.)

Dimensions:

25 in. x 17.2 in. x 17.3 in.
(63.5 cm. x 43.7 cm. x 43.9 cm.)

Total Unit Weight:

- Camera Only: 32.8 lbs. (14.1 kg.)
- 17-inch monitor - 39.5 lbs. (17.9 kg.)
- 20-inch monitor - 41 lbs. (18.6 kg.)
- 22-inch monitor - 42.1 lbs. (19.1 kg.)
- 24-inch monitor - 43.7 lbs. (19.8 kg.)

Controls:

- On / Off toggle switch
- Color mode select dial
- Magnification 16-position dial
- Freeze Frame push button
- Brightness dial
- Find push button
- Focus Lock and position locator push button
- Masks / Lines select push switch
- Masks / Lines position wheels
- Toggle PC and camera views push button

Power Supply:

Single wall plug, 50-60 Hz , 100-240V AC



Keep your place and make reading long texts easier with masking and reading lines

Designed & Assembled in the USA 



Freedom Scientific[®]

everything in sight.[™]

Call to Learn More
800-444-4443

www.FreedomScientific.com

Freedom Scientific is the world's leading manufacturer of assistive technology products for those who are low vision or blind.

Visit us online at www.FreedomScientific.com to learn more about these and other innovative products from Freedom Scientific that make living with low vision easier.



RUBY[®]
Pocket-Sized Handheld Video Magnifier for maximum portability



MAGic[®]
HD Screen Magnification Software with built-in screen reading features



ONYX[®] Deskset XL
Three-way Portable Video Magnifier for easy magnification near or far