

# MAGic Keystrokes

## General

Start MAGic	ALT+CTRL+M
Start MAGic Pro Scripting Edition	INSERT+O
Minimize/Restore the MAGic User Interface	MAGic Key+U
Pass Through Key	MAGic Key+3
Refresh Screen	MAGic Key+ESC
Open MAGic Context Menu	MAGic Key+APPLICATION Key

## Magnification, Views, & Enhancements

Increase Magnification	MAGic Key+NUM PAD PLUS or CTRL+EQUALS
Decrease Magnification	MAGic Key+NUM PAD MINUS or CTRL+DASH
Magnification Toggle On/Off	MAGic Key+DELETE
Toggle between 1x Magnification	MAGic Key+SHIFT+DELETE
Switch to Next Magnified View	MAGic Key+NUM PAD SLASH
Switch to Previous Magnified View	SHIFT+MAGic Key+NUM PAD SLASH
Toggle Color Enhancements On/Off	MAGic Key+F12
Toggle Cursor Enhancements On/Off	MAGic Key+F8
Toggle Mouse Enhancements On/Off	MAGic Key+F4

## Panning and Screen Movement

Pan Right	MAGic Key+RIGHT ARROW
Pan Left	MAGic Key+LEFT ARROW

Pan Down	<b>MAGic Key+DOWN ARROW</b>
Pan Up	<b>MAGic Key+UP ARROW</b>
Increase Speed during Panning	<b>NUM PAD PLUS</b>
Decrease Speed during Panning	<b>NUM PAD MINUS</b>
Move to Left Edge of Screen in Magnified View	<b>MAGic Key+HOME</b>
Move to Right Edge of Screen in Magnified View	<b>MAGic Key+END</b>
Move to top of Screen	<b>MAGic Key+PAGE UP</b>
Move to Bottom of Screen	<b>MAGic Key+PAGE DOWN</b>
Move to Beginning of Next Line	<b>MAGic Key+ENTER</b>

## Tracking

Tracking Toggle On/Off	<b>MAGic Key+R</b>
Move View to Mouse Pointer	<b>MAGic Key+NUM PAD 5</b>
Move View to Cursor	<b>SHIFT+MAGic Key+NUM PAD 5</b>
Tether Mouse to Display	<b>CTRL+INSERT+NUM PAD MINUS</b>

## Cursor Movement

Route Cursor to Mouse	<b>INSERT+NUM PAD PLUS</b>
Route Mouse to Cursor	<b>INSERT+NUM PAD MINUS</b>

## Locator Mode

Start Locator	<b>MAGic Key+F10</b>
Stop Locator and Move to Selected Area	<b>ENTER</b>
Stop Locator without Moving the Magnified Area	<b>ESC</b>

Move Left	LEFT ARROW
Move Right	RIGHT ARROW
Move Up	UP ARROW
Move Down	DOWN ARROW
Move to End of Line	END
Move to Beginning of Line	HOME
Move to Top of Screen	PAGE UP
Move to Bottom of Screen	PAGE DOWN

## Quick View Frame

Toggle Quick View Frame On/Off	MAGic Key+Q
Set upper left corner of frame	MAGic Key+LEFT BRACKET
Set lower right corner of frame	MAGic Key+RIGHT BRACKET
Create frame around current dialog box option	MAGic Key+LEFT BRACKET twice quickly
Open Quick View Manager	MAGic Key+9

## Speech

Silence Current Speech	CTRL
Enable/Disable Speech	MAGic Key+F3
Language Toggle	MAGic Key+F11
Set Screen Echo	INSERT+S
Select Verbosity Scheme	ALT+INSERT+S

# Document and Text Reading

Read All from Location	<b>MAGic Key+SPACEBAR or INSERT+DOWN ARROW</b>
Read Active Window (reads from the beginning of document window)	<b>MAGic Key+SPACEBAR (twice) or INSERT+DOWN ARROW (twice)</b>
Say Currently Selected Text	<b>MAGic Key+CTRL or INSERT+SHIFT+DOWN ARROW</b>
Say Current Character	<b>NUM PAD 5</b>
Say Current Word	<b>MAGic Key+ALT or INSERT+NUM PAD 5</b>
Spell Current Word	<b>MAGic Key+ALT (twice) or INSERT+NUM PAD 5 (twice)</b>
Say Next Word	<b>INSERT+RIGHT ARROW</b>
Say Prior Word	<b>INSERT+LEFT ARROW</b>
Select Next Word	<b>INSERT+SHIFT+RIGHT ARROW</b>
Select Prior Word	<b>INSERT+SHIFT+LEFT ARROW</b>
Say Current Line	<b>MAGic Key+L or INSERT+UP ARROW</b>
Say Next Line	<b>DOWN ARROW</b>
Say Prior Line	<b>UP ARROW</b>
Read to Cursor	<b>INSERT+HOME</b>
Read from Cursor	<b>INSERT+PAGE UP</b>
Say Current Sentence	<b>MAGic Key+S</b>
Say Next Sentence	<b>CTRL+MAGic Key+S</b>
Say Current Paragraph	<b>MAGic Key+P</b>
Read Column from Current Position	<b>MAGic Key+C</b>
Say Contents of Clipboard	<b>MAGic Key+X</b>

**Note:** When using the document and text reading commands, if you move the mouse **MAGic** reads from the current location of the mouse pointer. If you move the insertion point, **MAGic** reads from the location of the insertion point.

## Dialog Box Reading

Say Title Bar	<b>MAGic Key+T or INSERT+T</b>
Read Bottom of Window	<b>INSERT+PAGE DOWN</b>
Say Top Line of Window	<b>INSERT+END</b>
Say Default Button in Dialog Box	<b>INSERT+E</b>
Say Window Prompt and Text	<b>MAGic Key+TAB or INSERT+TAB</b>
Say Current Access Key	<b>SHIFT+NUM PAD 5</b>

## Table Reading

Say Current Cell	<b>ALT+CTRL+NUM PAD 5</b>
Say Cell to Right	<b>ALT+CTRL+RIGHT ARROW</b>
Say Cell to Left	<b>ALT+CTRL+LEFT ARROW</b>
Say Cell Below	<b>ALT+CTRL+DOWN ARROW</b>
Say Cell Above	<b>ALT+CTRL+UP ARROW</b>
Say First Cell	<b>ALT+CTRL+HOME</b>
Say Last Cell	<b>ALT+CTRL+END</b>
Say First Cell in Column	<b>ALT+CTRL+SHIFT+UP ARROW</b>
Say Last Cell in Column	<b>ALT+CTRL+SHIFT+DOWN ARROW</b>
Say First Cell in Row	<b>ALT+CTRL+SHIFT+LEFT ARROW</b>
Say Last Cell in Row	<b>ALT+CTRL+SHIFT+RIGHT ARROW</b>

## Typing Echo

Cycle through Indicate Capitalization options	<b>MAGic Key+F2</b>
Cycle through Typing Echo Unit options	<b>INSERT+2</b>

## Mouse Echo

Mouse Echo Toggle	<b>MAGic Key+M</b>
-------------------	--------------------

**Note:** Toggling the mouse echo rotates between speak lines, speak words, and mouse echo off.

Some actions in MAGic do not have hot keys assigned. You can define hot keys for these features from within the **Keyboard Manager** dialog box.