

TOPAZ[®] PHD

Portable Video Magnifier



Designed for Your Lifestyle!

Intuitive and fully-featured video magnifier provides brilliant image quality in a compact and portable design!



*folds easily
for transport*



Features

- ▶ Lightweight and foldable
- ▶ Easy-to-use, bright controls
- ▶ Customizable high-contrast colors
- ▶ Save to SD Card
- ▶ PC Connectivity

Learn more about getting back to the things that matter most.

www.FreedomScientific.com



1-866-569-8474 • 727-803-8000 (outside US)

TOPAZ[®] PHD

Portable Video Magnifier



Enjoy Freedom and Independence Again

Discover the freedom of independently reading and viewing the details that matter... in the home, office, classroom, or on the go.



SPECIFICATIONS*

MAGNIFICATION

12-in monitor: 1.7x to 24x

15-in monitor: 2.3x to 31x

DIMENSIONS H/W/D

13.5 in / 12.5 in / 2.4 in
34.3 cm / 31.7 cm / 6.1 cm

13.5 in / 16.0 in / 2.4 in
34.3 cm / 40.6 cm / 6.1 cm

WEIGHT

8.4 lbs / 3.8 kg w/o battery
9.4 lbs / 4.3 kg w/battery

9.4 lbs / 4.3 kg w/o battery
10.4 lbs / 4.7 kg w/battery

STANDARD COLOR MODES

True full color High-contrast full color
Black on white Yellow on blue
White on black Yellow on black
+27 Programmable

POWER SUPPLY

Input voltage 100-240V AC
Frequency 50-60 Hz
Single wall plug

CAMERA

14 MP Sensor

**subject to change*

WORKING HEIGHT UNDER CAMERA

9.75 in / 24.8 cm

STANDARD

PC connectivity Self-view camera
SD card slot Carrying case

OPTIONS

XY table
Lithium-ion battery pack



www.FreedomScientific.com

For product details, including advanced options, specifications, and warranty, please visit our website or scan this code.

750543-001 Rev. D



1-866-569-8474 • 727-803-8000 (outside US)