

***Full-featured video magnifier
in a compact, portable design***



The TOPAZ Ultra offers the same reading comfort, ease of use, and functionality as a desktop video magnifier - but in a compact and foldable design.



The ideal solution for individuals who need a lightweight magnification solution in different locations - at home, in the office or classroom, and on the go.



Use the TOPAZ Ultra right out of the box to magnify text and objects up close, in the distance, and even look at yourself, with its convenient 3-in-1 camera.

Enjoy crisp, clear, full-HD images on your choice of a 15.6- or 17.3-inch screen. The monitor height is fully adjustable, ensuring you always have the right reading angle. Our unique folding design makes it extremely easy to take the TOPAZ Ultra with you, or to store it when you do not need it.

Features

- Full 1080p HD with autofocus
- Large and intuitive control buttons
- Video magnifier with large 15.6- or 17.3-inch screen that easily folds away
- Read text and view objects up close, in the distance and view yourself
- Portable with long, 5-hour battery life
- Unique zoom and pan features
- 3 selectable magnification ranges to suit your eye condition
- Read full pages without moving reading material thanks to the movable camera
- View larger images by connecting to your TV screen via built-in HDMI port